


● ***Autoclavable Jar & Lid Assembly***



Canister Specifications

Volume (ml):	1200
Material:	Polycarbonate
Max. Diameter (mm):	177.2
Min. Diameter (mm):	131.6
Height (mm):	185
Measurement Tolerances:	±10%
Max. Vacuum:	0.1MPa (750 mmHg)
Filter Efficiency:	≥99.9%
Packaging:	48pcs/Carton; 0.15CBM

Autoclavable Jar and Lid Assembly

The Kaiya Medical Autoclavable Jar and Lid Assembly is manufactured with high density polypropylene. The jar and lid is used with the Kaiya Medical Models SU01A Suction Unit. This is a reusable jar and lid.