

## Kaiya 5L Oxygen Concentrator

Model - ZY-5AC

An ultra-quiet, FDA approved 5L oxygen concentrator with built in oxygen purity monitor sensor device, to measure the oxygen purity level at all time, it provide continuous 24 hour non-stop oxygen supply with purity up to 95% for home, hospital, clinical, nursing centre use.



Oxygen Purity Sensor



Digital Display Panel



FDA Approved



## Specifications

Power Supply	AC 115V / 60Hz or 220V / 50Hz
Power Consumption	< 340W or < 220 W
Output Pressure	40 ~ 60 Kpa
Output Flow	5 liter per minute
Oxygen Purity	93% ± 3%
Sound Level	< 50dB (A)
Classification	Class II Type B
Dimensions	36 cm x 29 cm x 58 cm
Weight	15.2 Kg / 12 Kg