

Compressor Nebulizer

Model - KYWH1004

A highly efficient, stable, easily maintained compressor nebulizer with extended life, able to deliver smaller particles of respiratory medicine for home, hospital, clinical and nursing center use.



Extended Life



Continuous Flow



MMAD 5 Microns or less



Specifications

Power Supply	AC 220V / 50Hz
Power Consumption	160VA
Max. Pressure	241 Kpa
Max. Flow	6 liter per minute
Nebulizing Pressure	75 Kpa
Sound Level	< 70dB (A)
Classification	Class II Type B
Dimensions	19.5 cm x 17.3 cm x 10.2 cm
Weight	1.86 Kg