

5 L Oxygen Concentrator Model - ZY-5AC Ultra Quiet

This ultra-quiet 5L FDA approved oxygen concentrator offers more peaceful nights and 24-hour service for the home-care patient



Real time O2 purity monitor



FDA approved



Ultra quite < 45dBA



Specifications

Power Supply	AC 115V/60Hz or 220V / 50Hz
Power Consumption	< 340 W
Output Pressure	40 ~ 60 Kpa
Max. Oxygen Flow	5 liter per minute
Oxygen Purity	93 % ± 3%
Sound Level	< 45 dB (A)
Classification	Class II Type B
Dimensions	36 cm x 29 cm x 58 cm
Weight	15.2 Kg